

DHRS3

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:1000 - 1:2000 IHC 1:50 - 1:100

Calculated MW:34kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DHRS3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SDR1; RDH17; Rsdr1; DD83.1; SDR16C1; retSDR1;

Background:

Short-chain dehydrogenases/reductases (SDRs), such as DHRS3, catalyze the oxidation/reduction of a wide range of substrates, including retinoids and steroids.

To place an order, please [Click HERE](#).

Catalog #:A7530

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9249

Isotype:IgG

Swiss Prot:O75911

Purity:Affinity purification

For research use only.