

DVL2

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:79kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human DVL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

t

Catalog #:A7551

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1856

Isotype:IgG

Swiss Prot:O14641

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the dishevelled (dsh) protein family. The vertebrate dsh proteins have approximately 40% amino acid sequence similarity with Drosophila dsh. This gene encodes a 90-kD protein that undergoes posttranslational phosphorylation to form a 95-kD cytoplasmic protein, which may play a role in the signal transduction pathway mediated by multiple Wnt proteins. The mechanisms of dishevelled function in Wnt signaling are likely to be conserved among metazoans.

To place an order, please [Click HERE](#).