

EIF4G1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 175kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human EIF4G1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

P220; EIF4F; EIF4G; EIF4GI; PARK18; EIF-4G1;

Catalog #: A7552

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1981

Isotype: IgG

Swiss Prot: Q04637

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a component of the multi-subunit protein complex EIF4F. This complex facilitates the recruitment of mRNA to the ribosome, which is a rate-limiting step during the initiation phase of protein synthesis. The recognition of the mRNA cap and the ATP-dependent unwinding of 5'-terminal secondary structure is catalyzed by factors in this complex. The subunit encoded by this gene is a large scaffolding protein that contains binding sites for other members of the EIF4F complex. A domain at its N-terminus can also interact with the poly(A)-binding protein, which may mediate the circularization of mRNA during translation. Alternative splicing results in multiple transcript variants, some of which are derived from alternative promoter usage.

To place an order, please [Click HERE](#).