

## INHBC

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:100 - 1:200

**Calculated MW:**38kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human INHBC

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

IHBC;

**Background:**

This gene encodes the beta C chain of inhibin, a member of the TGF-beta superfamily. This subunit forms heterodimers with beta A and beta B subunits. Inhibins and activins, also members of the TGF-beta superfamily, are hormones with opposing actions and are involved in hypothalamic, pituitary, and gonadal hormone secretion, as well as growth and differentiation of various cell types.

*To place an order, please [Click HERE](#).*

**Catalog #:**A7559

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**3626

**Isotype:**IgG

**Swiss Prot:**P55103

**Purity:**Affinity purification

For research use only.