

## PBX1

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 47kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human PBX1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

μg/ml

**Catalog #:** A7567

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 5087

**Isotype:** IgG

**Swiss Prot:** P40424

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a nuclear protein that belongs to the PBX homeobox family of transcriptional factors. Studies in mice suggest that this gene may be involved in the regulation of osteogenesis, and required for skeletal patterning and programming. A chromosomal translocation, t(1;19) involving this gene and TCF3/E2A gene, is associated with pre-B-cell acute lymphoblastic leukemia. The resulting fusion protein, in which the DNA binding domain of E2A is replaced by the DNA binding domain of this protein, transforms cells by constitutively activating transcription of genes regulated by the PBX protein family. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**To place an order, please [Click HERE](#).**