

ZSCAN26

Reactivity: Human Mouse

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:200

Calculated MW: 55kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ZSCAN26

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ZNF187; SRE-ZBP;

To place an order, please [Click HERE](#).

Catalog #: A7571

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7741

Isotype: IgG

Swiss Prot: Q16670

Purity: Affinity purification

For research use only.