

STK19

Reactivity:Human Mouse Rat

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:200

Calculated MW:41kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human STK19

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

G11; RP1; D6S60; D6S60E; HLA-RP1;

Background:

This gene encodes a serine/threonine kinase which localizes predominantly to the nucleus. Its specific function is unknown; it is possible that phosphorylation of this protein is involved in transcriptional regulation. This gene localizes to the major histocompatibility complex (MHC) class III region on chromosome 6 and expresses two transcript variants.

To place an order, please [Click HERE](#).

Catalog #:A7574

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8859

Isotype:IgG

Swiss Prot:P49842

Purity:Affinity purification

For research use only.