

NINL

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:156kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human NINL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NLP; dJ691N24.1;

To place an order, please [Click HERE](#).

Catalog #:A7584

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:22981

Isotype:IgG

Swiss Prot:Q9Y2I6

Purity:Affinity purification

For research use only.