

CACNG4

Reactivity:Human Mouse Rat

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:200

Calculated MW:37kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CACNG4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A7588

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:27092

Isotype:IgG

Swiss Prot:Q9UBN1

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a type I transmembrane AMPA receptor regulatory protein (TARP). TARPs regulate both trafficking and channel gating of the AMPA receptors. This gene is part of a functionally diverse eight-member protein subfamily of the PMP-22/EMP/MP20 family and is located in a cluster with two family members, a type II TARP and a calcium channel gamma subunit.

To place an order, please [Click HERE](#).