

DIMT1

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:35kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DIMT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DIM1; DIMT1L; HUSSY5; HSA9761;

To place an order, please [Click HERE](#).

Catalog #:A7589

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:27292

Isotype:IgG

Swiss Prot:Q9UNQ2

Purity:Affinity purification

For research use only.