

BANP

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 56kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BANP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BEND1; SMAR1; SMARBP1;

Catalog #: A7595

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54971

Isotype: IgG

Swiss Prot: Q8N9N5

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that binds to matrix attachment regions. The protein forms a complex with p53 and negatively regulates p53 transcription, and functions as a tumor suppressor and cell cycle regulator. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).