

TSSK3

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:30kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TSSK3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TSK3; SPOGA3; STK22C; STK22D;

Catalog #:A7605

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:81629

Isotype:IgG

Swiss Prot:Q96PN8

Purity:Affinity purification

For research use only.

Background:

This gene encodes a kinase expressed exclusively in the testis that is thought to play a role in either germ cell differentiation or mature sperm function.

To place an order, please [Click HERE](#).