

TTBK2

Reactivity:Human Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50

Calculated MW:137kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TTBK2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TTBK; SCA11;

Background:

This gene encodes a serine-threonine kinase that putatively phosphorylates tau and tubulin proteins. Mutations in this gene cause spinocerebellar ataxia type 11 (SCA11); a neurodegenerative disease characterized by progressive ataxia and atrophy of the cerebellum and brainstem.

To place an order, please [Click HERE](#).

Catalog #:A7609

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:146057

Isotype:IgG

Swiss Prot:Q6IQ55

Purity:Affinity purification

For research use only.