

STK32A

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:46kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human STK32A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

YANK1;

To place an order, please [Click HERE](#).

Catalog #:A7613

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:202374

Isotype:IgG

Swiss Prot:Q8WU08

Purity:Affinity purification

For research use only.