

CENPK

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 32kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CENPK

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P33; Solt; CENP-K; FKSG14; AF5alpha;

Catalog #: A7627

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64105

Isotype: IgG

Swiss Prot: Q9BS16

Purity: Affinity purification

For research use only.

Background:

CENPK is a subunit of a CENPH (MIM 605607)-CENPI (MIM 300065)-associated centromeric complex that targets CENPA (MIM 117139) to centromeres and is required for proper kinetochore function and mitotic progression.

To place an order, please [Click HERE](#).