

## CENPQ

---

**Reactivity:** Human Mouse

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 31kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human CENPQ

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CENP-Q; C6orf139;

**Catalog #:** A7628

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 55166

**Isotype:** IgG

**Swiss Prot:** Q7L2Z9

**Purity:** Affinity purification

For research use only.

**Background:**

CENPQ is a subunit of a CENPH (MIM 605607)-CENPI (MIM 300065)-associated centromeric complex that targets CENPA (MIM 117139) to centromeres and is required for proper kinetochore function and mitotic progression.

*To place an order, please [Click HERE](#).*