

UBIAD1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 37kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human UBIAD1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SCCD; TERE1;

Catalog #: A7640

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29914

Isotype: IgG

Swiss Prot: Q9Y5Z9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein thought to be involved in cholesterol and phospholipid metabolism.

Mutations in this gene are associated with Schnyder crystalline corneal dystrophy.

To place an order, please [Click HERE](#).