

## CLPTM1

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 76kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human CLPTM1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:** A7656

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 1209

**Isotype:** IgG

**Swiss Prot:** O96005

**Purity:** Affinity purification

For research use only.