

DLX6

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:20kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DLX6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A7667

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1750

Isotype:IgG

Swiss Prot:P56179

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of a homeobox transcription factor gene family similar to the Drosophila distal-less gene. This family is comprised of at least 6 different members that encode proteins with roles in forebrain and craniofacial development. This gene is in a tail-to-tail configuration with another member of the family on the long arm of chromosome 7.

To place an order, please [Click HERE](#).