

## YBX1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:100 - 1:200

**Calculated MW:**36kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human YBX1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

jk

**Synonym:**

YB1; BP-8; CSDB; DBPB; YB-1; CSDA2; NSEP1; NSEP-1; MDR-NF1;

*To place an order, please [Click HERE](#).*

**Catalog #:**A7704

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**4904

**Isotype:**IgG

**Swiss Prot:**P67809

**Purity:**Affinity purification

For research use only.