

## PLCG1

---

**Reactivity:**Human

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**149kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human PLCG1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

q

**Synonym:**

PLC1; NCKAP3; PLC-II; PLC148; PLCgamma1;

**Catalog #:**A7711

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5335

**Isotype:**IgG

**Swiss Prot:**P19174

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene catalyzes the formation of inositol 1,4,5-trisphosphate and diacylglycerol from phosphatidylinositol 4,5-bisphosphate. This reaction uses calcium as a cofactor and plays an important role in the intracellular transduction of receptor-mediated tyrosine kinase activators. For example, when activated by SRC, the encoded protein causes the Ras guanine nucleotide exchange factor RasGRP1 to translocate to the Golgi, where it activates Ras. Also, this protein has been shown to be a major substrate for heparin-binding growth factor 1 (acidic fibroblast growth factor)-activated tyrosine kinase. Two transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*