

PPP3CC

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:58kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PPP3CC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

i

Synonym:

CNA3; CALNA3; PP2Bgamma;

Catalog #:A7714

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5533

Isotype:IgG

Swiss Prot:P48454

Purity:Affinity purification

For research use only.

Background:

Calcineurin is a calcium-dependent, calmodulin-stimulated protein phosphatase involved in the downstream regulation of dopaminergic signal transduction. Calcineurin is composed of a regulatory subunit and a catalytic subunit. The protein encoded by this gene represents one of the regulatory subunits that has been found for calcineurin. Three transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).