

## MAP2K3

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**39kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human MAP2K3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

iors

**Synonym:**

MEK3; MKK3; MAPKK3; PRKMK3; SAPKK2; SAPKK-2;

**Catalog #:**A7718

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5606

**Isotype:**IgG

**Swiss Prot:**P46734

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a dual specificity protein kinase that belongs to the MAP kinase kinase family. This kinase is activated by mitogenic and environmental stress, and participates in the MAP kinase-mediated signaling cascade. It phosphorylates and thus activates MAPK14/p38-MAPK. This kinase can be activated by insulin, and is necessary for the expression of glucose transporter. Expression of RAS oncogene is found to result in the accumulation of the active form of this kinase, which thus leads to the constitutive activation of MAPK14, and confers oncogenic transformation of primary cells. The inhibition of this kinase is involved in the pathogenesis of Yersina pseudotuberculosis. Multiple alternatively spliced transcript variants that encode distinct isoforms have been reported for this gene.

*To place an order, please [Click HERE](#).*