

RAC1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 21kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RAC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1µg

Synonym:

MIG5; Rac-1; TC-25; p21-Rac1;

Catalog #: A7720

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5879

Isotype: IgG

Swiss Prot: P63000

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).