

CLIP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 162kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CLIP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RSN; CLIP; CYLN1; CLIP170; CLIP-170;

Catalog #: A7722

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6249

Isotype: IgG

Swiss Prot: P30622

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene links endocytic vesicles to microtubules. This gene is highly expressed in Reed-Sternberg cells of Hodgkin disease. Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).