

SLC6A1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:67kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human SLC6A1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MAE; GAT1; GABATR; GABATHG;

Background:

The SLC6A1 gene encodes a gamma-aminobutyric acid (GABA) transporter, which removes GABA from the synaptic cleft (Hirunsatit et al., 2009)

To place an order, please [Click HERE](#).

Catalog #:A7730

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6529

Isotype:IgG

Swiss Prot:P30531

Purity:Affinity purification

For research use only.