

## TFAP2C

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**49kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human TFAP2C

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ERF1; TFAP2G; hAP-2g; AP2-GAMMA;

**Background:**

The protein encoded by this gene is a sequence-specific DNA-binding transcription factor involved in the activation of several developmental genes. The encoded protein can act as either a homodimer or heterodimer with other family members and is induced during retinoic acid-mediated differentiation. It plays a role in the development of the eyes, face, body wall, limbs, and neural tube.

*To place an order, please [Click HERE](#).*

**Catalog #:**A7739

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**7022

**Isotype:**IgG

**Swiss Prot:**Q92754

**Purity:**Affinity purification

For research use only.