

TRPS1

Reactivity:Human Mouse

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW:142kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TRPS1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GC79; LGCR;

Catalog #:A7743

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7227

Isotype:IgG

Swiss Prot:Q9UHF7

Purity:Affinity purification

For research use only.

Background:

This gene encodes a transcription factor that represses GATA-regulated genes and binds to a dynein light chain protein. Binding of the encoded protein to the dynein light chain protein affects binding to GATA consensus sequences and suppresses its transcriptional activity. Defects in this gene are a cause of tricho-rhino-phalangeal syndrome (TRPS) types I-III.

To place an order, please [Click HERE](#).