

WNT7B

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:39kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human WNT7B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A7746

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7477

Isotype:IgG

Swiss Prot:P56706

Purity:Affinity purification

For research use only.

Background:

This gene is a member of the WNT gene family, which consists of structurally related genes that encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. Among members of the human WNT family, this gene product is most similar to WNT7A protein.

To place an order, please [Click HERE](#).