

ST8SIA2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:42kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ST8SIA2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

STX; SIAT8B; HsT19690; ST8SIA-II;

Catalog #:A7748

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8128

Isotype:IgG

Swiss Prot:Q92186

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a type II membrane protein that is thought to catalyze the transfer of sialic acid from CMP-sialic acid to N-linked oligosaccharides and glycoproteins. The encoded protein may be found in the Golgi apparatus and may be involved in the production of polysialic acid, a modulator of the adhesive properties of neural cell adhesion molecule (NCAM1). This protein is a member of glycosyltransferase family 29.

To place an order, please [Click HERE](#).