

DGKE

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:64kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DGKE

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DGK; DAGK5; DAGK6; NPHS7;

Background:

Diacylglycerol kinases are thought to be involved mainly in the regeneration of phosphatidylinositol (PI) from diacylglycerol in the PI-cycle during cell signal transduction. When expressed in mammalian cells, DGK-epsilon shows specificity for arachidonyl-containing diacylglycerol. DGK-epsilon is expressed predominantly in testis.

To place an order, please [Click HERE](#).

Catalog #:A7752

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8526

Isotype:IgG

Swiss Prot:P52429

Purity:Affinity purification

For research use only.