

NAT8

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 26kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NAT8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GLA; CML1; CCNAT; Hcml1; ATase2; TSC501; TSC510;

Catalog #: A7759

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9027

Isotype: IgG

Swiss Prot: Q9UHE5

Purity: Affinity purification

For research use only.

Background:

This gene, isolated using the differential display method to detect tissue-specific genes, is specifically expressed in kidney and liver. The encoded protein shows amino acid sequence similarity to N-acetyltransferases. A similar protein in *Xenopus* affects cell adhesion and gastrulation movements, and may be localized in the secretory pathway. A highly similar paralog is found in a cluster with this gene.

To place an order, please [Click HERE](#).