

PDE4DIP

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:265kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PDE4DIP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MMGL; CMYA2;

Background:

The protein encoded by this gene serves to anchor phosphodiesterase 4D to the Golgi/centrosome region of the cell. Defects in this gene may be a cause of myeloproliferative disorder (MBD) associated with eosinophilia. Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).

Catalog #:A7765

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9659

Isotype:IgG

Swiss Prot:Q5VU43

Purity:Affinity purification

For research use only.