

ARPC3

Reactivity:Human Mouse Rat

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW:21kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ARPC3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARC21; p21-Arc;

Catalog #:A7767

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10094

Isotype:IgG

Swiss Prot:O15145

Purity:Affinity purification

For research use only.

Background:

This gene encodes one of seven subunits of the human Arp2/3 protein complex. The Arp2/3 protein complex has been conserved through evolution and is implicated in the control of actin polymerization in cells. Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).