

RRAGA

Reactivity: Human Mouse Rat

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW: 37kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RRAGA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FIP1; RAGA; FIP-1;

To place an order, please [Click HERE](#).

Catalog #: A7771

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10670

Isotype: IgG

Swiss Prot: Q7L523

Purity: Affinity purification

For research use only.