

## KDM5B

---

**Reactivity:** Human

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 176kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human KDM5B

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CT31; PLU1; PUT1; PLU-1; JARID1B; PPP1R98; RBBP2H1A;

*To place an order, please [Click HERE](#).*

**Catalog #:** A7772

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 10765

**Isotype:** IgG

**Swiss Prot:** Q9UGL1

**Purity:** Affinity purification

For research use only.