

CFHR3

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:37kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CFHR3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FHR3; HLF4; CFHL3; FHR-3; DOWN16;

Background:

The protein encoded by this gene is a secreted protein, which belongs to the complement factor H-related protein family. It binds to heparin, and may be involved in complement regulation. Mutations in this gene are associated with decreased risk of age-related macular degeneration, and with an increased risk of atypical hemolytic-uremic syndrome. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).

Catalog #:A7775

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10878

Isotype:IgG

Swiss Prot:Q02985

Purity:Affinity purification

For research use only.