

STARD3

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:51kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human STARD3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CAB1; es64; MLN64;

Background:

This gene encodes a member of a subfamily of lipid trafficking proteins that are characterized by a C-terminal steroidogenic acute regulatory domain and an N-terminal metastatic lymph node 64 domain. The encoded protein localizes to the membranes of late endosomes and may be involved in exporting cholesterol. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).

Catalog #:A7776

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10948

Isotype:IgG

Swiss Prot:Q14849

Purity:Affinity purification

For research use only.