

CKAP4

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:66kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human CKAP4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

p63; CLIMP-63; ERGIC-63;

To place an order, please [Click HERE](#).

Catalog #:A7777

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10970

Isotype:IgG

Swiss Prot:Q07065

Purity:Affinity purification

For research use only.