

FAF1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 74kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human FAF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

hFAF1; CGI-03; HFAF1s; UBXD12; UBXN3A;

Catalog #: A7779

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11124

Isotype: IgG

Swiss Prot: Q9UNN5

Purity: Affinity purification

For research use only.

Background:

Interaction of Fas ligand (TNFSF6) with the FAS antigen (TNFRSF6) mediates programmed cell death, also called apoptosis, in a number of organ systems. The protein encoded by this gene binds to FAS antigen and can initiate apoptosis or enhance apoptosis initiated through FAS antigen. Initiation of apoptosis by the protein encoded by this gene requires a ubiquitin-like domain but not the FAS-binding domain.

To place an order, please [Click HERE](#).