

PACSIN2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 56kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PACSIN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SDPII;

Catalog #: A7781

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11252

Isotype: IgG

Swiss Prot: Q9UNF0

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the protein kinase C and casein kinase substrate in neurons family. The encoded protein is involved in linking the actin cytoskeleton with vesicle formation by regulating tubulin polymerization. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).