

SUZ12

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:83kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SUZ12

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

CHET9; JJAZ1;

Catalog #:A7786

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23512

Isotype:IgG

Swiss Prot:Q15022

Purity:Affinity purification

For research use only.

Background:

This zinc finger gene has been identified at the breakpoints of a recurrent chromosomal translocation reported in endometrial stromal sarcoma. Recombination of these breakpoints results in the fusion of this gene and JAZF1. The protein encoded by this gene contains a zinc finger domain in the C terminus of the coding region.

To place an order, please [Click HERE](#).