

COG4

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:89kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human COG4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

COD1; CDG2J;

Background:

The protein encoded by this gene is a component of an oligomeric protein complex involved in the structure and function of the Golgi apparatus. Defects in this gene may be a cause of congenital disorder of glycosylation type IIj. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).

Catalog #:A7792

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:25839

Isotype:IgG

Swiss Prot:Q9H9E3

Purity:Affinity purification

For research use only.