

SAMHD1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 72kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SAMHD1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DCIP; CHBL2; HDDC1; MOP-5; SBB188;

Catalog #: A7794

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 25939

Isotype: IgG

Swiss Prot: Q9Y3Z3

Purity: Affinity purification

For research use only.

Background:

This gene may play a role in regulation of the innate immune response. The encoded protein is upregulated in response to viral infection and may be involved in mediation of tumor necrosis factor-alpha proinflammatory responses. Mutations in this gene have been associated with Aicardi-Goutieres syndrome.

To place an order, please [Click HERE](#).