

EHF

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:35kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human EHF

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ESE3; ESEJ; ESE3B;

Catalog #:A7795

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:26298

Isotype:IgG

Swiss Prot:Q9NZC4

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein that belongs to an ETS transcription factor subfamily characterized by epithelial-specific expression (ESEs). The encoded protein acts as a transcriptional repressor and may be involved in epithelial differentiation and carcinogenesis. Three transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).