

MLH3

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 164kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MLH3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HNPCC7;

Catalog #: A7796

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27030

Isotype: IgG

Swiss Prot: Q9UHC1

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the MutL-homolog (MLH) family of DNA mismatch repair (MMR) genes. MLH genes are implicated in maintaining genomic integrity during DNA replication and after meiotic recombination. The protein encoded by this gene functions as a heterodimer with other family members. Somatic mutations in this gene frequently occur in tumors exhibiting microsatellite instability, and germline mutations have been linked to hereditary nonpolyposis colorectal cancer type 7 (HNPCC7). Several alternatively spliced transcript variants have been identified, but the full-length nature of only two transcript variants has been determined.

To place an order, please [Click HERE](#).