

DKK4

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DKK4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DKK-4;

Catalog #: A7797

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27121

Isotype: IgG

Swiss Prot: Q9UBT3

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that is a member of the dickkopf family. The secreted protein contains two cysteine rich regions and is involved in embryonic development through its interactions with the Wnt signaling pathway. Activity of this protein is modulated by binding to the Wnt co-receptor and the co-factor kremen 2.

To place an order, please [Click HERE](#).