

SCG3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 53kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SCG3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

SGIII;

Catalog #: A7799

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29106

Isotype: IgG

Swiss Prot: Q8WXD2

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. Granins may serve as precursors for biologically active peptides. Some granins have been shown to function as helper proteins in sorting and proteolytic processing of prohormones; however, the function of this protein is unknown. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).