

TNNI3K

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 93kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TNNI3K

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CARK; CCDD;

Catalog #: A7802

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51086

Isotype: IgG

Swiss Prot: Q59H18

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that belongs to the MAP kinase kinase kinase (MAPKKK) family of protein kinases. The protein contains ankyrin repeat, protein kinase and serine-rich domains and is thought to play a role in cardiac physiology.

To place an order, please [Click HERE](#).