

## SNX7

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 45kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human SNX7

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:** A7805

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 51375

**Isotype:** IgG

**Swiss Prot:** Q9UNH6

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region like some family members, and its exact function is unknown. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 11.

*To place an order, please [Click HERE](#).*